



Scapa Soft-Pro[®] Polyurethane Gel

Redefining Skin Friendly



Scapa
Healthcare

Scapa Soft-Pro[®] Polyurethane Gel

Scapa Soft-Pro[®] Polyurethane Gel combines the ultimate in breathability, security and skin-friendly comfort and is ideal for a variety of skin contact applications.

When combined with a high MVTR backing the breathable properties of Scapa Soft-Pro[®] Polyurethane Gel allows for moisture to pass away from the skin, reducing the risk of maceration.

Skin Friendly Redefined: Secure, Low Trauma, Breathable

- High MVTR
 - Reduce risk of maceration
- Minimal trauma upon removal
- Minimal trauma during wear
- Waterproof
- Suitable for Gamma Sterilization





Image represents Scapa Soft-Pro[®] 6052 Polyurethane Gel with Polyurethane Carrier.

Secure, breathable and skin friendly adhesive


Scapa Soft-Pro[®] Polyurethane Gel provides a breathable and skin friendly adhesive when a secure fixation is required on active patients.

MVTR Breathability*

g/m²/24 hour – upright cup
(results of adhesive coated on PU film)



Peel Strength*



Scapa Soft-Pro[®] Polyurethane Gel

Product	Carrier	Adhesive	Liner	Used in:
6052	25 micron PU Film with PE Carrier	100 g/m ² PU Gel	White release paper	Medical Devices, Consumer Wellness, Advanced Wound Care

* Test data on file.


Why Scapa?

- We care about quality
- Our high-performing materials will differentiate your products
- We're experts in skin friendly adhesives
- We can manage everything, from product design to delivery
- We're trusted worldwide

Our expertise in product development means we are the ideal business partner in the healthcare industry. Our goal is to provide you with a reliable, trusted, regulated, safe and cost-effective partnership.

Skin Friendly

Turn-Key Solutions



Think.
Create.
Deliver.

scapahealthcare.com

sales.healthcare@scapa.com

Scapa North America

Tel: +1 860 688 8000

Scapa Europe

Tel: +44 (0)161 301 7400

Scapa Asia

Tel: +603 5569 2529

For technical and application enquiries please email: r&d@scapa.com

© 2014 Scapa Group plc. Scapa Soft-Pro® is a registered trademark of Scapa Group plc