



Scapa Soft-Pro® Silicone Gel

Apply the gentle touch



Scapa
Healthcare

Scapa Soft-Pro® Silicone Gel


Scapa Soft-Pro® Silicone Gel is highly conformable and easy to apply and remove. The silicone gel layer ensures secure fixation and gentle removal from skin. Scapa Soft-Pro® Silicone Gel is an ideal fixation solution for patients with fragile and sensitive skin. It is suitable for patients requiring repeated application and removal over the same area.

Scapa Soft-Pro® Silicone Gel may be incorporated into your final product design, providing skin friendly adhesion for the patient. Scapa Soft-Pro® Silicone Gel can also be used for fixation of medical devices and can be produced in double-sided and single-sided variants.

Ideal for use on fragile and sensitive skin:

- Low trauma, minimizes pain on removal
- Use on open or intact skin
- Repositionable
- Conforms to body contours

Adhesive comparison chart





Scapa Soft-Pro® Silicone is available in standard validated adhesion strengths or may be customized for specific applications

| Product | Width (mm / inch) | Coat Weight (g/m ²) | Liner | Film | Adhesive Colour | Silicone Gel | | Acrylic Adhesive | |
|---------|-------------------|---------------------------------|-------|------|-----------------|----------------------------------|-----------------------------|--------------------------------|--------------------------------|
| | | | | | | Securement Strength [#] | Liner Release (N/25mm, 90°) | Peel (N/25mm, Stainless Steel) | Liner Release (N/25mm, T-Peel) |
| RX1423S | 635 / 25 | 150 | PC | TPU | Clear | Medium | 0.2 | 10 | 0.15 |
| 6054 | 1120 / 57 | 150 | LDPE | PU | Clear | High | 0.40 | 9 | 0.17 |
| 6051 | 1120 / 57 | 200 | LDPE | PU | Clear | High | 0.2 | 13 | 0.2 |

Perforating

Scapa has the capability to perforate Silicone Gel. Perforation of the silicone gel allows for transmission of fluid and exudate minimizing the risk of skin maceration. It is ideal for use with absorbable components in a layered construction.

Perforation dimensions




Why Scapa?

- We care about quality
- Our high-performing materials will differentiate your products
- We're experts in skin friendly adhesives
- We can manage everything, from product design to delivery
- We're trusted worldwide

Our expertise in product development means we are the ideal business partner in the healthcare industry. Our goal is to provide you with a reliable, trusted, regulated, safe and cost-effective partnership.

Skin Friendly

Turn-Key Solutions



Think.
Create.
Deliver.

scapahealthcare.com
sales.healthcare@scapa.com

Scapa North America
Tel: +1 860 688 8000

Scapa Europe
Tel: +44 (0)161 301 7400

Scapa Asia
Tel: +603 5569 2529

For technical and application enquiries please email: r&d@scapa.com

© 2014 Scapa Group plc. Scapa Soft-Pro® is a registered trademark of Scapa Group plc