

RenPro™ 188

Professional Trainers Tape



Strength in the heat of battle...



RENFREW
ATHLETICS

1.800.542.9715 • renfrewathletics@scapa.com • www.renfrewathletics.com

PROFESSIONAL TRAINERS TAPE

RenPro™ 188

RenPro™ 188 Professional Trainers Tape is a tight-weave premium bleached cotton cloth coated with a zinc oxide-impregnated rubber adhesive. This full-featured athletic trainers tape is specially formulated to meet the demands of professional athletes and deliver top performance under intense conditions.

- Bleached cotton cloth
- Strong zinc oxide-impregnated rubber adhesive
- Resists failure caused by moisture
- Consistent unwind
- Easy tear for speed-wrapping



8.5 X 7.625

Renfrew 106

Foam Pre-wrap



Unbeatable Performance...



RENFREW
ATHLETICS

1.800.542.9715 • renfrewathletics@scapa.com • www.renfrewathletics.com

FOAM PRE-WRAP

Renfrew 106

Renfrew 106 Foam Pre-wrap is a lightweight foam under-wrap that provides additional comfort for the athlete and facilitates easy removal of all Renfrew Athletics trainers tapes.

- “Open-cell” polyurethane foam provides breathability
- Absorbent foam keeps moisture from outer wrap



Available in eight colors:


Black, Red, Purple, Orange, Tan, Blue, Yellow and Pink.

8.5 X 8



Renfrew 287

PolyCotton Trainers Tape



Sturdy Performance...

RENFREW
ATHLETICS

1.800.542.9715 • renfrewathletics@scapa.com • www.renfrewathletics.com

POLYCOTTON TAPE

Renfrew 287

Renfrew 287 is an unbleached polyester/cotton blended cloth, coated with a zinc oxide impregnated rubber adhesive. This durable and cost effective alternative to bleached tapes delivers all the performance and characteristics required of a professional trainers tape without the added costs involved with bleaching.

- Natural unbleached polyester / cotton cloth
- Strong zinc oxide-impregnated rubber adhesive
- Resists failure caused by moisture
- Fully-functional cost-effective alternative to bleached tapes



8.5 X 8.375

RenFrew 136

Wrestling Mat Tape



Instant Tape Down...

RENFREW
ATHLETICS

1.800.542.9715 • renfrewathletics@scapa.com • www.renfrewathletics.com

WRESTLING MAT TAPE

RenFrew 136

Renfrew 136 Wrestling Mat Tape is a heavy-duty polyethylene tape used for seaming and securing wrestling mats. Its tough backing resists tears and withstands the stresses of intense use over extended periods, and its rubber adhesive bonds extremely well to mat surfaces. This product is very conformable and its excellent tear characteristics make it easy to apply and ideal for securing the outer edges of mats without curling.



8.5 X 9.125

RenKon™ 208

Heavyweight Elastic Stretch Tape



208

One Tough Competitor...

RENFREW
ATHLETICS

1.800.542.9715 • renfrewathletics@scapa.com • www.renfrewathletics.com

HEAVYWEIGHT ELASTIC STRETCH TAPE

RenKon™ 208

RenKon™ Heavyweight Elastic Stretch Tape is a narrow-woven cotton cloth, coated with a porous zinc oxide-impregnated rubber adhesive that provides excellent adhesion and easy unwind. The soft edges of the narrow-weave cloth provide increased comfort for the athlete. RenKon is ideal for injury prevention and compression management.

- 100% high-twist cotton yarns for maximum conformability and consistent stretch
- Porous adhesive provides breathable properties
- Fast and easy unwind
- Narrow-weave cloth provides soft edges and increased comfort
- Tan colored with red centerline



Scapa

8.5 X 9.875

RenLite™ 209

Lightweight Stretch Tape



Tight to the Finish...

**RENFREW
ATHLETICS**

1.800.542.9715 • renfrewathletics@scapa.com • www.renfrewathletics.com

LIGHTWEIGHT STRETCH TAPE

RenLite™ 209

RenLite™ 209 Lightweight Stretch Tape is a highly conformable cotton cloth, coated with a porous zinc oxide-impregnated rubber adhesive that provides secure adhesion and easy unwind. The product's high tensile strength offers maximum support and makes it an ideal product for injury prevention and directional support.

- Highly conformable cotton cloth
- Porous adhesive provides breathable properties
- High tensile strength provides maximum support
- Fast and easy unwind
- Easy tear and no-slip grip for high-speed taping



8.5 X 10.25



Renfrew PFT™ 210

Patella-femoral Tape



Tough in the Trenches...

RENFREW
ATHLETICS

1.800.542.9715 • renfrewathletics@scapa.com • www.renfrewathletics.com

PATELLA-FEMORAL TAPE

Renfrew PFT™ 210

Renfrew PFT™ 210 Patella-femoral Tape is a heavy-duty cotton/rayon cloth, coated with an aggressive zinc oxide-impregnated rubber adhesive. The high-tensile rigid cloth provides maximum support and limits extreme ranges of motion.

- Strong rigid cloth for maximum support
- Aggressive adhesive ensures tape stays in place for the duration of activity
- High-tensile cloth limits motion
- Pinked edges for easy tearing and high-speed taping



8.5 X 10.625

